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Silikon-Foley-Katheter

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Silicone Foley Type Catheter

Description

All catheters basically comprise a 2-lumen or 3-lumen shaft with proximal funnel, inflation valve and distal retaining balloon with various tip designs. Balloon fill volumes are in ml. and the shaft sizes are Charrière (Ch.) and are indicated on the funnel of each individual catheter and on the individual pack label and outer carton. The distal tip type is indicated on the funnel of each individual pack label and/or outer carton.

The MPI-balloon catheter is made from 100% silicone.

Size range

From 6 Ch. to 26 Ch.

Length will vary between 20 and 45 cm.

Indications

These catheters are indicated for routine drainage of the bladder or for routine post-operative drainage and irrigation of the bladder.

Read all instructions carefully before use.

Preparation for use

Remove catheter from its sterile pouch using aseptic technique. Visually inspect the catheter for any mechanical damage. If a trial inflation of the balloon is carried out, ensure the balloon is empty prior to insertion. Only sterile water should be used to fill the balloon.

Directions for use

Lubricate the catheter using only a water based lubricant or gel (with or without local anesthetic as required). Using sterile technique as indicated in published literature; insert the catheter fully into the urethra ensuring that the balloon is beyond the bladder neck. All angled catheters tips correspond to the direction of the inflation valve. A Luer syringe may be used to inflate the balloon. Inflate the balloon with sterile water. With some products a pre-filled syringe is enclosed. The volume of the balloon is indicated on the funnel of the catheter or/and on the label of the packaging.

Warning: Do not over inflate.

Ensure function by checking for urine flow from the drainage funnel. To deflate the balloon prior to removal, insert a Luer syringe as appropriate. Aspirate the syringe gently to remove the inflation fluid. DO NOT use excessive aspiration on the syringe during deflation as this may cause a vacuum collapse of the inflation lumen, which may impair aspiration of inflation fluid. Maximum recommended indwelling time is 42 days.

The inflation balloon holds the catheter in the bladder. If the balloon loses fluid, e.g. by diffusion of the filling fluid through the balloon membrane, the catheter can, with time, slide out of the bladder accidentally. In order to prevent this danger, the inflation liquid should be withdrawn from the catheter at least once every 2 weeks and the balloon should then be immediately reinflated with the stated nominal filling volume. The patient should be routinely monitored in accordance with accepted procedures and the catheter removed after a suitable interval, not longer than 42 days, as

determined by a physician or other suitably qualified personnel.

Potential complications

If difficulty is encountered aspirating the balloon with a syringe, a rare and infrequent reported event, the branch of the catheter with the valve should be cut with sharp scissors at the bifurcation or the balloon ruptured according to established procedures reported in medical literature. Should it be necessary to rupture the balloon, care must be taken to remove all fragments from the patient. Incorrectly positioned catheters can cause urethral damage if the balloon is inflated within the urethra.

Irritation of the urethral mucosa, blockage of the catheter due to encrustation and catheter-induced infections are documented complications with some patients.

Caution

For urological use only. Federal law restricts this device to sale by or on the order of a physician (USA & Canada).

Use only a water based lubricant.

Always inflate the balloon sterile water with glycerin 9:1 / USP 29 purified water.

Never clamp the catheter; a catheter plug or valve may be used if required.

Never use a syringe to puncture the catheter shaft for urine sampling.

There is a danger of needle stick injury and the catheter may be compromised.

Presentation

Each product is supplied in an individual sterile pouch.

The product is guaranteed sterile unless the individual packaging has been opened or damaged.

This product is for single use only. Do not reuse or resterilize. Re-use may cause a risk of infection and/or compromise functional reliability.

Single use product

These products are for single use only. Reprocessing may result in degraded performance or loss of functionality. Validated cleaning and sterilization methods and instructions for reprocessing to original specifications are not available for these devices. This product is not designed to be cleaned, disinfected or sterilized.

Storage

Do not store our products in extreme temperatures or wet conditions.

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