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
## Content

### Matthys Drain®

- Catheter puncture kit
  - Perforation unit (A) (perforating cannula with protective cap, catheter with plug and protective sheath), Ø 2.7 mm, length 450 mm
  - Slip-on fixture (B) (slip-on cannula with protective cap, 3-way stopcock, one-way valve)
  - Anti-kink sleeve (C)
  - Flexible Y non-return valve unit (D)
  - Syringe 30 ml
  - Secretion bag 2000 ml

Materials used: PE, PVC, ABS, PC, PP, stainless steel, PUR.

### Fields of application:

Drainage of transpleural exudate and exudate of the pleural and pulmonary cavities; abscesses, pericardial effusion, ascites, treatment of pneumothorax, irrigation of pleural empyema and other body cavities, drainage of subcutaneous emphysema, instillation of medication into body cavities. Transtracheal administration of oxygen in emergencies.

**Literature:** please request as desired

**Contraindications:** coagulation defects (absolute or relative contraindication)

### Risks:

The known risks exist according to the type of application; in particular bleeding, injury to the visceral pleura and lungs, damage to the lungs through excessive suction (only suction equipment with reliable suction restriction should be used), circulatory disturbances on extraction of excessive amounts of fluid in pleural effusion and ascites; infections, pericardial effusion, pneumothorax through incorrect application. The patient should be informed of the risks!!

### Warning!

- Do not use if the packaging is damaged.
- Do not resterilise (disposable product)
- Must be replaced for each patient. Otherwise, there is a risk of cross-contamination which causes severe problems for the patient.
- Do not use in children, during pregnancy or while breastfeeding.

### Caution!


The catheter puncture kit should be used only under sterile conditions once the sterile packaging has been opened.

The secretion bag should be hung or laid below the level of the patient, for only in this way is a slow drainage of secretion possible without suction equipment.


Duration of application: The catheter material can remain in situ for a maximum of 30 days. If there are signs of local or general infection which may have been caused by the catheter, it should be removed and a fresh catheter laid.

### Directions for use for the physician

- Handling the kit and laying the catheter must take place under sterile conditions after a thorough disinfection of the hands. The indicated point of perforation should be anaesthetised locally.
- Remove perforation unit from the packaging.
- Remove the protective cap from the perforation unit.
- With careful twisting movements pierce the skin with the perforation cannula and catheter and insert it into the body cavity.


- Push the catheter in its protective sheath through the perforation cannula into the body cavity to a level of not less than the last calibration ring. Hold the perforation cannula with one hand and with the other hand grasp the catheter through the protective sheath and push in as far as required.
- Hold the catheter with its protective sheath and withdraw the perforation cannula from the skin perforation.

### Caution!

On removal of the perforation cannula the catheter must not be drawn out of the body any further than its first black (wide) mark. Withdraw the catheter as required, observing the marks (black marks at intervals from 12 to 22 mm). Never pull the catheter out through the perforation cannula to avoid cutting it off.

- When the tip of the perforation cannula is approx. 2 cm outside the body, the catheter must be held with two fingers or fixed with a plaster strip so that the perforation cannula and protective sheath can be completely removed via the catheter.
- Remove the slip-on fixture from the packaging and check that the 3-way stopcock is closed to the slip-on cannula (B).


- Slip the anti-kink sleeve onto the end of the catheter that is sealed with the plug.
- Flex the catheter approx. 2 cm from the plug, so that it is sealed, and hold it with two fingers.
- Remove the plug from the end of the catheter and connect the 3-way stopcock of the slip-on cannula in its place. Push in as far as possible.

### Caution!

Make sure that the catheter is firmly connected to the slip-on cannula (push together until there is an overlap of at least 1 cm)

- Push the anti-kink sleeve back over the catheter end and onto the slip-on cannula.
- Sample material can be removed or medication introduced via the side connection of the 3-way stopcock, using the 30 ml syringe.
- In the case of a pneumothorax the use of the flexible Y non-return valve unit (D) is particularly recommended. This should be attached with the long length of tube to the non-return valve of the connecting unit.
- Firmly screw the secretion bag (only in the case of effusion) on to the non-return valve of the flexible Y non-return valve unit.
- Set the 3-way stopcock so that the connection is open. Check with the 30 ml syringe at the short length of tube of the Y non-return valve unit (F) whether air or fluid can be pumped off. The two non-return valves prevent aspiration and/or the backflow of any substance from/into the body cavity.

- If there is spontaneous drainage of fluid, the secretion bag is placed or hung below the point of perforation. In such a case there is usually no need for pumping off with the 30 ml syringe or for connecting any suction equipment, (max. vacuum  $9 \times 10^{-4}$  Pa).
- If air and fluid are being drained, it is recommended that the secretion bag be perforated so that the fluid collects in the bag and the air can escape (e.g. in the case of a seropneumothorax).

Negative pressure occurs in the pleural space during puncturing of the pleural cavity with the catheter perforation unit, and a small pneu-

mothorax ensues corresponding to the amount of air in the airtight catheter sheath. This enables the pleural viscera to move out of the way of the needle point and also minimizes the risk of injury from the slightly protruding catheter tip. If there is excess pressure in the pleural space or a pleural effusion, the catheter sheath fills with air (tension pneumothorax) and/or secretion. In all of these cases the catheter must be pushed forward through the cannula at least as far as the first (wide) black graduation mark so that no air from outside can enter the pleural cavity through holes at the side, or fluid escape.

When installing a continuous suction appliance we recommend attaching the vacuum pump to the Luer-Lock connection of the flexible Y non-return valve unit. The catheter should be covered with sterile slitted gauze and either stitched to the body or attached with fixing plaster.

### Caution!

The position of the drain should be verified by means of an X-ray. In a pneumothorax the patency of the entire system should be checked by emerging the Y non-return valve unit (connection E) in a glass of water whilst the patient coughs. If not patent (no escape of air bubbles), the system is flushed with sterile saline solution injected through the syringe connection (F) of the 3-way stopcock in the direction of the patient and/or secretion bag.

### Medication Instillation

- Close 3-way stopcock (B) to catheter.
- Remove cap from free port of 3-way stopcock.
- Connect syringe with medication to 3-way stopcock.
- Open 3-way stopcock to connect syringe to catheter.
- Instilling medication.
- Close the three-way stopcock towards the catheter before removing the syringe.
- Remove the syringe and close with cap.

### Storage:

- Room temperature 15 - 40°C
- Relative humidity 18 - 80%
- Keep away from direct light!
- With appropriate storage of the product (avoiding hard or angular objects), use by the date shown on the packaging.

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Mediland GmbH