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Uretersplint directable

Description

Uretersplints consist of polyurethane, have a closed end or central opening for coaxial application. They are partially perforated throughout the entire length and marked every cm. X-ray. Splints are connected to the pusher and can be disconnected inside the bladder by turning the lock mechanism.

Length between the curves: 22 / 24 / 26 / 28 cm

Ch.	4,8	6,0	7,0	8,0
Max. diameter of guide wire [inch]	0,032	0,035	0,038	0,038

Indications

Directable splints are intended for use as routine ureter stents. Application has to be done by a physician or suitably qualified medical personnel.

Contraindications

Insurmountable ureteric obstruction.

Preparation of use

Length and diameter must be individually suitable to the patient. It is important that the outer diameter of the splint is compatible to the inner diameter of the endoscope. Remove the splints from their sterile pouches using an aseptic technique. Do not use kinked guide wires.

Application

A. Axial application (splint with guide wire)

1. The stent is connected to the pusher and the guide wire is partly premounted. It needs to be pushed to the very tip of the stent.
2. Fix the stent on the pusher with a clamp.
3. Introduce the assembled stent, renal end first, into the working channel of the in-situ cystoscope or ureterscope and gently place into the ureter under x-ray control.
4. Remove the clamp from the pusher as required to allow further progress of the stent or pusher into the endoscope channel.
5. Once the placement of the renal tip has been confirmed in the renal pelvis using established techniques, the guide wire should be slowly withdrawn to allow the renal coil to form.
6. Now dislock the pusher by separating the LL-connection by turning. After dislocking remove the introducer.
7. The removal of the stent may be done by a grasping forceps through a cystoscope locating the vesical end of the splint.
8. The stents may also be placed antegrade via percutaneous nephrostomy.

B. Co-axial application (first introducing the guide wire before pushing the splint)

1. In case of coaxial application push the splint carefully over the in-situ guide wire. Then push the splint with the pre-mounted introducer over the guide wire until it has its final position.
2. When the location of the splint tip in the renal

- pelvis has been verified by established techniques then, slowly withdraw the guide wire. The tip will have the correct coil form then.
3. Now dislock the introducer by turning the LL-connection.
4. The removal of the stent may be done by a grasping forceps through a cystoscope locating the vesical end of the splint.
5. The splints can also be placed antegrade through percutane nephrostomy.

Potential complications

Infections by improper use, blockade of the stent due to encrustation, perforation and irritation of bladder and ureter.

In order to assure the flow in the stent the patient should be routinely monitored.

Supply

Every set is supplied sterile in an individual pouch.

The sterility of the set and its components is guaranteed as long as the packaging unit is unopened and intact, but no longer than stated on the label.

Storage

The product should be stored in a dry place and protected against direct an indirect light, moisture and heat sources.

Single-use product

These products are for single use only. Reprocessing may result in degraded performance or loss of functionality. Validated cleaning and sterilization methods and instructions for reprocessing to original specifications are not available for these devices. This product is not designed to be cleaned, disinfected or sterilized.

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