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CE 1254

## EN

### Insertion stylet for tracheal tubes

### Instructions for disinfection, cleaning and sterilization

#### Disinfection and cleaning:

Factory-new instruments should – as far as they are not sterile packed – be washed prior to sterilization. Previously used instruments must be cleaned, even after having been disinfected. Thermic disinfection should be performed at an average temperature of 93°C.

For the cleaning of our instruments we recommend using cleaning agents that are suitable for medical instruments consisting of rubber, latex and plastics. Such cleaning agents are available from dealers specialized in medical technology products. We point out that the manufacturers instructions are to be followed. Hard brushes and other materials that can damage the surface of the instruments should not, under any circumstance, be used.

Hot air treatment is harmful due to free flow of oxygen to the material surface, and can cause premature aging. The drying temperature should, therefore, not exceed 70 °C (158 °F). It is recommended that the dried articles be removed from the hot air section as soon as possible.

When using washing machines featuring temperature restriction in their drying program, care should be taken not to exceed 70 °C (158 °F). The drying temperature should be checked from time to time.

During treatment with water or aqueous solvents, small amounts of water can enter the materials of many of our instruments. In this case, the instruments occasionally show a milky white haze, which can be removed by drying at a higher temperature. Drying can take from 30 minutes up to 2 hours, whereby the temperature may not exceed 70 °C (158 °F). Upon storage in a suitable environment the milky white haze will disappear by itself. Softening of thermoplastic instruments attributable to heating will disappear again upon cooling.

Neodisher and similarly aggressive cleaning agents will dissolve softeners out of PVC and will result in increased rigidity of the plastic coating. This is unpreventable, however, must be respected for the lifetime of the instrument ( visible by hardening of the end sections ).

#### Steam sterilization:

121°C, 20 minutes without vacuum drying  
134°C, 5 minutes without vacuum drying

Attention should be paid to the fact that each steam sterilization will accelerate natural aging.

#### Gas sterilization:

Our instruments can, in principle, be subjected to sterilization with ethylene oxide gas or similar gases, or gas mixtures.

Please ensure that the ethylene oxide residual concentration should be less than 4 mg/product prior to approval for use.

#### General instructions of the manufacturer:

The insertion stylets are manufactured from a tension-free extruded PVC tube, which is thermally sealed at the tip. The flexible metal core consists of stainless steel 1.4104 and is secured by an open PVC splint to allow unimpeded pressure compensation during reprocessing.

The insertion stylets were factory-tested on 10 application cycles and subsequently sterilized according to the following procedure:

1. Mechanical alkaline cleaning and disinfection at 93°C with subsequent steam sterilization at 121°C during 20 minutes.
2. Mechanical alkaline cleaning and disinfection at 93°C with subsequent steam sterilization at 134°C during 5 minutes.

#### Important:

- During conditioning, the instruments must be positioned upright!!
- No automatic vacuum drying!!

**Automatic vacuum drying following the sterilization process can damage the PVC coating. Water that has entered the PVC coating under the influence of steam will abruptly evaporate under vacuum and will destroy the hot and softened PVC layer.**

Our instruments can incur damages that can lead to complete destruction when coming into contact with the following compounds:

- All fats and oils, such as Vaseline and paraffin oil, lubricants with or without anaesthetising action
- Organic solvents, such as petroleum ether, benzene, ethers, esters, ketones, chlorinated hydrocarbons, etc.
- Oxidising compounds such as chlorine or peracetic acid should not be used in an acidic pH range.
- Disinfecting agents that contain phenol or phenolic compounds.

**We recommend maximum 5 conditionings - instruments that, after frequent use, have become markedly discoloured, show cracking or have developed a sticky surface, exhibit hardening or other changes, are defective and are not to be used for safety reasons.**

These instructions on disinfection, cleaning and sterilization are provided to the best of our knowledge. We, however, decline any liability derived from this.

MEDICOPLAST INTERNATIONAL GmbH  
Heusweilerstrasse 100  
66557 Illingen  
Germany

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